

UCLA 2P Miniscope Assembly Guide

Component Inventory and Validations

The first step in assembly is to validate all components are sourced and ready to be used. For the headpiece, it's straightforward to think about breaking the parts down into their respective subassemblies and ensure that all parts are assembled. Please find links to each component either within the GitHub repository or respective manufacturer. Some components, like the electronic PCBs, the custom lenses made by Optics Technology, and diced filters, are custom made or modified and need to be either constructed yourself or with the help of manufacturers. We are hoping to put together batch orders for components such that users can receive all needed parts for fabrication in a single shipment. For more information please contact pgolshani@mednet.ucla.edu

These components are required for the **Main Housing Assembly**:

Component	Part Number	Manufacturer	Unit	Price	Subtotal	Link
Main Housing	Main_Housing.STL	Custom	1	\$1.00	\$1.00	Main Housing
Lens Cover	Lens_Cover.STL	Custom	1	\$1.00	\$1.00	Lens Cover
Flex-PCB	FLEX-PCB	Custom	1	\$600.00	\$600.00	Flex-PCB
Scan Lens	AC050-010-B	ThorLabs	2	\$51.38	\$102.76	Scan Lens
Tube Lens*	OT-1 (Custom)	Optics Technology	1	\$1,100.00	\$1,100.00	Tube Lens
2P Dichroic (Diced)	ZT775sp-2p, 8x8.5x1mm	Chroma	1	\$540.00	\$540.00	2P Dichroic
Torx Screws	96817A704	McMaster-Carr	4	\$0.63	\$2.53	Torx Screws
Optical Glue	NOA68	Norland	1	\$40.70	\$40.70	Optical Glue
					\$2,387.99	
* Tube lens is custom fabricated, please contact pgolshani@mednet.ucla.edu for more information						

The following components are required for the **Detector Assembly**:

Component	Part Number	Manufacturer	Unit	Price	Subtotal	Link
Detector Housing	Detector_Housing.STL	Custom	1	\$1.00	\$1.00	Detector Housing
Emission Filter Support	Emission_Filter_Support.STL	Custom	1	\$1.00	\$1.00	Emission Support
Collection Lens	47-895	Edmund Optics	1	\$62.50	\$62.50	Collection Lens
Emission Filter (Cut)*	ET750sp, 4mm OD x1mm	Chroma	2	\$150.00	\$300.00	Emission Filter
1P Dichroic (Diced)**	T550lpxr, 4x4x1mm	Chroma	1	\$112.50	\$112.50	1P Dichroic
SiPM Detector	S13360-3075PE	Hamamatsu	2	\$52.88	\$105.76	SiPM Detector
					\$582.76	
* The emission filter used is a modified form of the ET750sp-2p8 (linked) to reduce the thickness						
** The cost per filter was reduced in our order due to multiple 4mm filters being produced from one standard filter after dicing						

The following components are required for the **Collimator Assembly**:

Component	Part Number	Manufacturer	Unit	Price	Subtotal	Link
Collimator Housing	Collimator_Housing.STL	Custom	1	\$1.00	\$1.00	Collimator Housing
Collimation Asphere	KGA170-B	Newport	1	\$107.00	\$107.00	Collimation Asphere
ETL	A-25H1	Varioptic	1	\$135.00	\$135.00	ETL
Torx Screws	96817A704	McMaster-Carr	4	\$0.63	\$2.53	Torx Screws
					\$245.53	

The following components are required for the **Lower Housing Assembly**:

Component	Part Number	Manufacturer	Unit	Price	Subtotal	Link
Lower Housing	Lower_Housing.STL	Custom	1	1	1	Lower Housing
Objective Lens*	OT-2 (Custom)	Optics Technology	1	4000	4000	Objective Lens
O-Ring		McMaster-Carr	1	0.1596	0.1596	O-Ring
					\$4,001.16	
* Tube lens is custom fabricated, please contact pgolshani@mednet.ucla.edu for more information						

The necessary **tools** for assembling the microscope are:

Component	Part Number	Manufacturer	Unit	Price	Subtotal	Link
UV-Curing System*	CS20K2	ThorLabs	1	\$2,634.66	\$2,634.66	UV-Curing System
Flathead Driver	84705260150	Wiha	1	\$5.17	\$5.17	Flathead Driver
Torx Screwdriver	52995A24	McMaster-Carr	1	\$9.48	\$9.48	Torx Driver
Scan Lens Reamer	8851A18	McMaster-Carr	1	\$24.18	\$24.18	Scan Lens Reamer
Tube Lens Reamer	8851A21	McMaster-Carr	1	\$28.78	\$28.78	Tube Lens Reamer
Ferrule Reamer	2777A128	McMaster-Carr	1	\$116.97	\$116.97	Ferrule Reamer
					\$2,819.24	
*This is the replacement for what we use in the lab (now obsolete part no. LED-200, ThorLabs)						

Please note that there are also custom interface circuits that are needed to build the complete system. They are covered in another document and are not referenced in this document further since the main purpose is to describe the method for constructing the microscope headpiece itself. Additionally, you will need an appropriate length of cleaved and sealed HC-920 optical fiber. For now, we will assume that the reader has appropriately launched the 2P laser to the HC-920 optical fiber, with sufficient efficiency.